


TT-COVER K-LIGHT™ MINI

TRAUMA



- WARMING
- HYGIENIC
- TIMESAVING



TT-Cover K-Light™ Mini is a disposable cover with separate head cushions with velcro. It is designed for optimum tensile strength. TT-Cover K-Light™ Mini is for paediatric use, 1 - 8 years.



ESON COMFORT
 Idé-Möbler i Landeryd AB
 Centrumvägen 7
 SE-314 94 LANDERYD

Tel +46 (0)371-40051
info@esoncomfort.com

Copyright 2016 Eson Comfort / Idé-Möbler i Landeryd AB.

All rights reserved

PRODUCT INFORMATION

TT-Cover K-Light™ Mini is a disposable cover with separate head cushions with velcro. It is developed to achieve durability and safe immobilisation. The fabrics is composed to reduce the loss of patientwarmth and to prevent contamination of the TraumaTransfer™-mattress and the trolley. It is designed for optimum tensile strength. Tested for a work load of 150 kg in wet condition. Destruction of the used product shall be done by burning to reinvest the energy contents.

NOTE! Contains no polyurethane that produces potential cancerogenic wast when burning.

ART.NO	NAME
TTCOVERKL	TT-Cover K-Light™ disposable cover for TraumaTransfer™ (25 pcs)
TTCOVERKLM	TT-Cover K-Light™ Mini disposable cover for TraumaTransfer™ Mini (25 pcs)
TTHEADFIXLH	Head cushion for TT-Cover K-Light™ (20 pcs) (order separately)
TTHEADFIXMLH	Head cushion for TT-Cover K-Light™ Mini (20 pcs) (order separately)

Extras / Spareparts

Special size, special shape.

CE Registrationcode 344-2757-118173

ESON COMFORT
www.esoncomfort.com
"WE INVENT THE FUTURE"