


ACCESSORIES 2020

TRAUMATRANSFER™


TRAUMA MATTRESS




1. TT4NG - TT4 Next Gen™ - Standard Body strap set for TraumaTransfer™. 4 pcs. of identical straps. For TraumaTransfer™ Mini it's the article TTM3NG and contains 3 identical straps. No exponents Velcro straps and padding of HD-plastics, serves for good hygiene.



2. TT44NG - TT44 Next Gen™ - Set with 4-point Body strap and 3 standard straps. For TraumaTransfer™ Mini it's the article TTM33NG and contains the 4-point Body strap and 2 standard straps. Approved for transport with TraumaTransfer™.



3. TTXNG - TTX Next Gen™ - Set with Cross sling and 3 standard straps. For TraumaTransfer™ Mini it's TTMXNG and contains the cross sling and 2 standard straps. Approved for transport with TraumaTransfer™.



4. BEDSLIDASET - Set with sliding boards
8 pcs. BedSlide™ Original, 4 pcs. BedSlide™ Large and 2 pcs. BedSlide™ Extra Large. Each size can also be ordered separately.



5. PELFIX - Immobilisation
Set with 2+2 polyethene cushions. **NOTE!** ALWAYS use immobilisation cushions together with TraumaTransfer™.



6. EASYCOLLV - AEROresc® - EASY Collar combines the best from the market leaders in collars. EASY Collar can be adjusted to 16 sizes, to obtain full adaption. Comes in an adult, military and paediatric version. It's also certified for, and a great complement to TraumaTransfer™.



7. TTTRANS - TT-Transporter™ levelling mattress
Used when patients are transported on TraumaTransfer™. Developed for use with Alfa-stretcher, with TT-EasySafety™ fixation.



8. TTEASYS - TT-EasySafety™ Safety solution. Set with 1 smaller and 1 wider belt solution. **NOTE!** Shall ALWAYS be used when patient is transported on TraumaTransfer™.



ESON COMFORT
Idé-Möbler i Landeryd AB
Centrumvägen 7
SE-314 94 LANDERYD

Tel +46 (0)371-40051
info@esoncomfort.com

Copyright 2016 Eson Comfort / Idé-Möbler i Landeryd AB.

All rights reserved



ESON COMFORT

“WE INVENT THE FUTURE”

www.esoncomfort.com