

ACCESSORIES TRAUMATRANSFER™


TRAUMA MATTRESS



1. **TT4 - Standard Body strap set** for TraumaTransfer™. 4 pcs. of identical straps. For TraumaTransfer™ Mini it's the article **TTM3** and contains 3 identical straps.



2. **TT44 - Set with 4-point Body strap** and 3 standard straps. For TraumaTransfer™ Mini it's the article **TTM33** and contains the 4-point Body strap and 2 standard straps. Approved for transport with TraumaTransfer™.



3. **TTX - Set with Cross sling** and 3 standard straps. For TraumaTransfer™ Mini it's **TTMX** and contains the cross sling and 2 standard straps. Approved for transport with TraumaTransfer™.




4. **BEDSLIDASET - Set with sliding boards** 8 pcs. **BedSlide™ Original**, 4 pcs. **BedSlide™ Large** and 2 pcs. **BedSlide™ Extra Large**. Each size can also be ordered separately.



5. **PELFIX - Immobilisation** Set with 2+2 polyethene cushions. **NOTE!** ALWAYS use immobilisation cushions together with TraumaTransfer™.



6. **SIP - Immobilisation** Set with 2+2+2 cushions. **NOTE!** Alternative to PELFIX with a special invented **AutoAdapt™** function. Follows the shape of the body with continuous immobilisation.



7. **TTTRANS - TT-Transporter™ levelling mattress** Used when patients are transported on TraumaTransfer™. Developed for use with Alfa-stretcher, with **TT-EasySafety™** fixation.



8. **TTEASYS - TT-EasySafety™** Safety solution. Set with 1 smaller and 1 wider belt solution. **NOTE!** Shall ALWAYS be used when patient is transported on TraumaTransfer™.



ESON COMFORT
Idé-Möbler i Landeryd AB
Centrumvägen 7
SE-314 94 LANDERYD

Tel +46 (0)371-40051
info@esoncomfort.com

Copyright 2016 Eson Comfort / Idé-Möbler i Landeryd AB.

All rights reserved



ESON COMFORT
"WE INVENT THE FUTURE"
www.esoncomfort.com