


TT-COVER K-LIGHT 2020™

TRAUMA



- WARMING
- HYGIENIC
- TIMESAVING



TT-Cover K-Light 2020™ is a disposable cover with separate head cushions. It is designed for optimum tensile strength. The picture to the right shows TT-Headfix K-Light 2020™.



ESON COMFORT
Idé-Möbler i Landeryd AB
Centrumvägen 7
SE-314 94 LANDERYD

Tel +46 (0)371-40051
info@esoncomfort.com

Copyright 2016 Eson Comfort / Idé-Möbler i Landeryd AB.

All rights reserved

PRODUCT INFORMATION

TT-Cover K-Light 2020™ is a disposable cover with separate head cushions. It is developed to achieve durability and safe immobilisation. The fabrics is composed to reduce the loss of patientwarmth and to prevent contamination of the TraumaTransfer™-mattress and the trolley. It is designed for optimum tensile strength. Tested for a work load of 150 kg in wet condition. Destruction of the used product shall be done by burning to reinvest the energy contents.

NOTE! Contains no polyurethane that produces potential cancerogenic wast when burning.

ART.NO	NAME
TTCOVERKL2	TT-Cover K-Light 2020™ disposable cover for TraumaTransfer™ (25 pcs)
TTCOVERKLM2	TT-Cover K-Light 2020™ Mini disposable cov for TraumaTransfer™ Mini (25 pcs)
TTHEADFIXKLH2	Head cushion for TT-Cover K-Light 2020™ (20 pcs) (order separately)
TTHEADFIXKLMH2	Head cushion for TT-Cover K-Light 2020™ (20 pcs) (order separately)

Extras / Spareparts

Special size, special shape.

CE Registrationcode 344-2757-118173

ESON COMFORT
www.esoncomfort.com
"WE INVENT THE FUTURE"