

TT-HEADFIX 2014™

TRAUMA



- SAFE
- DURABLE
- EASY TO USE



TT-Headfix 2014™ is a disposable head support for TraumaTransfer™. It is designed for optimum stability. It complements TT-Cover 2014™, when the whole cover isn't needed.



ESON COMFORT
Idé-Möbler i Landeryd AB
Centrumvägen 7
SE-314 94 LANDERYD

Tel +46 (0)371-40051
info@esoncomfort.com

Copyright 2016 Eson Comfort / Idé-Möbler i Landeryd AB.

All rights reserved

PRODUCT INFORMATION

TT-Headfix™ is a disposable head support. It is developed to achieve durability and safe immobilisation. The same solution is an important part of the TT-Cover 2014™ concept. It is designed for optimum stability. Tested for an angle load of 60 kg. Destruction of the used product shall be done by burning to reinvest the energy contents.

NOTE! Contains no polyurethane that produces potential cancerogenic waste when burning.

ART.NO
TTHEADFIX2

NAME
TT-Headfix 2014™ head-immobilisation for TraumaTransfer™ (20pcs.)

Extras / Spareparts

Special size, special shape.

Swedish patent search: 0402834-6
Publication: 529895IPC-klass: A61G 1/04, A61F 5/37
CE Registrationcode 344-2757-75049UNSPSC Code: 42.17.17.01GMDN Code: 35181

ESON COMFORT
www.esoncomfort.com
"WE INVENT THE FUTURE"