

TT-HEADFIX K-LIGHT™

TRAUMA



- SAFE
- DURABLE
- EASY TO USE



TT-Headfix K-Light™ is a disposable head support with separate head cushions, for TraumaTransfer™. It is designed for optimum stability. It complements TT-Cover K-Light™, when the whole cover isn't needed.



ESON COMFORT
 Idé-Möbler i Landeryd AB
 Centrumvägen 7
 SE-314 94 LANDERYD

Tel +46 (0)371-40051
info@esoncomfort.com

Copyright 2016 Eson Comfort / Idé-Möbler i Landeryd AB.

All rights reserved

PRODUCT INFORMATION

TT-Headfix K-Light™ is a disposable head support with separate head cushions. It is developed to achieve durability and safe immobilisation. The same solution is an important part of the TT-Cover K-Light™ concept. It is designed for optimum stability. Destruction of the used product shall be done by burning to reinvest the energy contents.

NOTE! Contains no polyurethane that produces potential cancerogenic waste when burning.

ART.NO

TTHEADFIXKL
 TTHEADFIXKLM
 TTHEADFIXKLLH
 TTHEADFIXKLMH

NAME

TT-Headfix K™ head-immobilisation for TraumaTransfer™ (20pcs.)
 TT-Headfix K™ Mini paediatric head-immobilisation for TT (20pcs.)
 Head cushion for TT-Headfix K-Light™ (20 pcs) (order separately)
 Head cushion for TT-Headfix K-Light™ Mini (20 pcs) (order separately)

Extras / Spareparts

Special size, special shape.

CE Registrationcode 344-2757-118173

ESON COMFORT
www.example.com

“WE INVENT THE FUTURE”