

TT-TRANSPORTER™

TRAUMA



- EASY TO USE
- SAFE
- DURABLE



TT-Transporter™ is a compensating mattress. It transform the angles in the chassis of Alfa-stretcher to a flat surface, during transport with TraumaTransfer™. Needs to be complemented with TT-EasySafety™ safety package to secure TraumaTransfer™ on top of TT-Transporter™ in the ambulance.



ESON COMFORT
Idé-Möbler i Landeryd AB
Centrumvägen 7
SE-314 94 LANDERYD

Tel +46 (0)371-40051
info@esoncomfort.com

Copyright 2016 Eson Comfort / Idé-Möbler i Landeryd AB.

All rights reserved

PRODUCT INFORMATION

The mattress is shaped to make a flat surface for TraumaTransfer™ on the Alfa-stretcher. The core has a volume weight of 44kg/m³. The coated fabrics is developed to be resistant to urine and exposure of cleaning-fluids with alcohol. Comes in a bag that protects the mattress and can store extras. Compatible with cigarett- and ignition test according to Swedish standard, without flame retardents. TT-Transporter™ comes with a 1 year technical warranty of function.

NOTE! Always use TT-Transporter™ together with TT-Easy Safety™.

ART.NO	NAME
TTTRAN	TT-Transporter mattress Black, Alfa-stretcher 49 x 192 x 10cm
TTSafe	TT-EasySafety safety package to TT-Transporter™/TraumaTransfer™
TT44	4-point safetybelt to TraumaTransfer™ (ALL pre-hospital transportation of patients on TraumaTransfer™ need to be secured with TT-Easy Safety™)

Extras/accessories

Specials

Cleaned with soap/alcohol

CE Registrationcode

Cleaning

ESON COMFORT
www.esoncomfort.com
"WE INVENT THE FUTURE"